RPE-IgG Conjugation Kit

Table of Contents

ntroduction	1
Package Contents and Storage Conditions	1
-eatures	2
Products and Contents	2
Additional Reagents Required	. 2
Shelf Life	2
gG Amount and Concentration and Buffers	. 2
RPE:IgG Molar Ratio	.3
Experimental Protocols	.3
RECOMMENDED ACCESSORIES	. 5

Introduction

R-Phycoerythrin is widely used as a fluorescent label in immunochemistry assays such as ELISA and in more complex techniques such as flow cytometry. Preparing bright, stable and reproducible antibody-PE conjugates is one of the biggest challenges of developing bead-based immunoassays and high quality reagents for flow cytometry. The OriGene RPE-IgG conjugation kit utilizes a novel chemistry to generate bright and highly reproducible RPE-IgG conjugates with a simple procedure. The resulting conjugates have been shown to be extremely stable, retaining 95% activity after storage for 30 days at 37°C with concentrations as low as $0.5 \mu g/mL$.

Package Contents and Storage Conditions

Kit Component	Storage Temp	Storage Notes
Concentrated Activator	-20°C	Keep the vial in the desiccated container as supplied in the kit
RPE-Z [™]	-20°C or 2-8°C	Does not need to be kept desiccated.
Quenching Reagent	-20°C or 2-8°C	Does not need to be kept desiccated.
Zeba Desalting Columnwith Collection Tube	nnwith 2-8°C Do not store at freezing	

Features

- Liquid-based reagents
- Completely scale-able: conjugate anywhere from 0.1 to 1 gram IgG per reaction.
- Supplies sufficient activated RPE to conjugate all IgG at a 1:1 RPE:IgG molar ratio.
- Highly efficient RPE incorporation purification not usually necessary.
- Customize the RPE:IgG ratio to create optimized conjugates for different applications.
- Conjugates have greatly improved stability vs Lightning-Link[™] and other chemistries.

Catalog Number	AR100080	AR100081	AR100082	AR100083	AR100084
For Labeling:	0.2 mg lgG	3x 100ug IgG	5 mg lgG	10 mg IgG	100 mg lgG
Concentrated Activator	10 µL	10 µL	10 µL	10 µL	20 µL
RPE-Z™ (16 mg/mL)	16 μL	36 µL	400 μL	0.80 mL	8.0 mL
Zeba Desalting Column	0.32 mg	3 each	8 mg	16 mg	160 mg
Quenching Reagent	25 μL	40 μL	60 µL	120 μL	1000 μL

Products and Contents

Additional Reagents Required But Not Supplied

- 1X Phosphate Buffered Saline (1X PBS)
- pH 7.2-7.5 Deionized water (dH2O)
- ml microcentrifuge tubes

Shelf Life

The performance of the product is guaranteed for a minimum of 12 months when stored as directed.

IgG Requirements

The IgG to be labeled should be at a minimum concentration of 0.8 mg/ml in pure 1X PBS and should not contain any preservatives or carriers such as sodium azide, Proclin 300 or BSA.

RPE:IgG Molar Ratio

This kit utilizes a 1:1 RPE:IgG molar ratio which is optimal for most conjugations reaction. However, lower or higher ratios may give better results depending upon the antibody characteristics and the intended end-use. To change the RPE:IgG molar ratio, vary the volume of RPE-Z[™] added to the conjugation reaction.

Conjugation Procedure - Overview





Before Beginning The Procedure

Remove the Concentrated Activator from the freezer. <u>Important: Allow sufficient time to let the</u> <u>container and contents come to room temperature before openingthe outer and inner vials.</u>

Note: The jar containing the Activator can be removed from the freezer up to 24hours before use.

Detailed Conjugation Procedure

- **1.** Measure the absorbance of the IgG solution at 280 nm using PBS as a blank.Divide the A280 by 1.40 to obtain the IgG concentration in mg/ml.
- 2. Dilute IgG to 1.20 mg/ml in 1X PBS (0.80 1.4 mg/ml is acceptable).
- 3. Add 100 μL of IgG solution to a new microcentrifuge tube.
- 4. Prepare a working dilution (1X) of Activator from Concentrated Activator indeionized water:

a. Add 2.0 uL of Concentrated Activator to 1300 µL of deionized water.

b. Immediately vortex to mix the solution thoroughly.

Note: The <u>1X</u> Activator must be used within 5 minutes of preparation. If more than5 minutes passes before use, discard the 1X Activator and prepare a fresh solution.

5. Add 2.0 uL of <u>1X</u> Activator to the 10 μ L aliquot of IgG and then mix thoroughlyby gentle vortexing.

6. Incubate the solution at room temperature for 1 hour.

Note: A longer incubation is not harmful and overnight incubations will besuccessful.

7. Desalt the complete 12 µL reaction volume into pure 1X PBS using the includedZeba spin column. See the attached desalting protocol.

Note: The activated IgG is stable and can be stored at 2-8°C for at least 4 months.

8. Add 2 µL of RPE-Z[™] to the desalted, activated IgG and mix by gentle vortexing.

9. Incubate the solution at room temperature for 2-24 hours.

ORIGENE

Note: Usable conjugates are produced after only 2 hours of conjugation. Larger and more potent conjugates will be produced after longer incubations.

10. Add 2 μL of Quenching Reagent to the reaction and mix by gentle vortexing.

11. Incubate the solution at room temperature for 1 hour.

Note: A longer incubation is not harmful and overnight incubations will besuccessful.

12. Test conjugate in the desired application.

Note: To improve conjugate performance, it may help to purify the conjugate from the unincorporated RPE and reaction components by size exclusion chromatography.

RECOMMENDED ACCESSORIES

For desalting IgG after activation - Order from ThermoFisher :

Sample Size	Description	Cat #
$2-12 \ \mu L$	Zeba Spin Desalting Columns, Micro (75µL), 7K MWCO	89877, 89878
30 - 130 μL	Zeba Spin Desalting Columns, 0.5 mL, 7K MWCO	89882, 89883
$200-700\;\mu L$	Zeba Spin Desalting Columns, 2 mL, 7K MWCO	89889, 89890
$500-2000 \ \mu L$	Zeba Spin Desalting Columns, 5 mL, 7K MWCO	89891, 89892
$700-4000 \ \mu L$	Zeba Spin Desalting Columns, 10 mL, 7K MWCO	89893, 89894

For concentrating IgG before or after activation or for concentrating the final conjugate – Order from MilliporeSigma:

Sample Size	Description	Cat #
Up to 500 µL	Amicon Ultra-0.5 Centrifugal Filter Unit with Ultracel-50 membrane	Z740176
Up to 2 mL	Amicon Ultra-2 Centrifugal Filter Unit with Ultracel-50 membrane	UFC205024
Up to 4 mL	Amicon Ultra-4 Centrifugal Filter Unit with Ultracel-50 membrane	UFC805008
Up to 15 mL	Amicon Ultra-15 Centrifugal Filter Unit with Ultracel-50 membrane	Z648000