



# Accel 3-in-1 EDTA Solution (20X Concentrate)

Storage: 2-8°C

Catalog No:	B22C-1L	1,000ml 🗌
-	B22C-500	500ml 🗌
	B22C-125	125ml
	B22C-50	50ml

## Description:

**Accel 3-in-1 EDTA Solution** is a revolutionary 3-in-1 pretreatment and antigen retrieval system intended for use in immunohistochemistry. This reagent is designed to deparaffinize, rehydrate, and unmask (HIER) formalin-fixed, paraffin-embedded tissue sections in one simple step prior to IHC staining using GBI Lab's Polymer Detection Systems. *This reagent can also be used as a conventional antigen retrieval buffer following traditional deparaffinization/hydration methods.* 

#### Kit Format:

liquid form (20X Concentrate).

Catalog Number	Description	Size
B22C-50	Accel 3-in-1 EDTA Solution (20X Concentrate)	50mL
B22C-125	Accel 3-in-1 EDTA Solution (20X Concentrate)	125mL
B22C-500	Accel 3-in-1 EDTA Solution (20X Concentrate)	500mL
B22C-1L	Accel 3-in-1 EDTA Solution (20X Concentrate)	1000mL

#### **Materials Required:**

- 1. Electric Pressure Cooker (we recommend the Cuisinart CPC-600 1000-Watt 6-Quart Electric Pressure Cooker) or Laboratory grade water bath
- 2. One metal rack
- 3. Two IHC staining dishes
- 4. Distilled water
- 5. PBS-Tween20 wash buffer

## **Electric Pressure Cooker Method:**

- 1. Prepare Ready-to-use (RTU) solution by diluting the concentrate 1:20 in distilled water.
- 2. Place up to 24 formalin-fixed, paraffin-embedded tissue slides in a slide rack.
- 3. Place slide rack into the first staining dish and fill with approximately 200 ml of RTU solution. (ensure tissue slides are fully immersed in RTU solution). Prepare a second staining dish filled with approximately 200 ml of RTU solution. The second staining dish will used later on as a rinse solution. *Staining dish size may vary, ensure sufficient quantity of RTU solution to cover the tissue slides.*
- 4. Place 1000 ml of distilled water in the base of pressure cooker unit and place a metal rack into the water. Make sure the bottom of the staining dish is not in direct contact with the bottom of the pressure cooker.
- 5. Place both staining dishes inside the pressure cooker. Lock pressure cooker lid in place. Set the pressure timer to 10-15 minutes. Make sure the vent switch is in the closed position. *The optimal incubation time should be determined by the user.*
- 6. Timer will start to count down when the correct pressure and temperature is reached.
- 7. After the timer finishes, turn the vent switch to release the pressure.
- 8. When all the pressure is released and you may safely remove the lid.
- 9. Remove both staining dishes from the pressure cooker. Wear proper safety equipment when preforming this step.
- 10. Agitate the slide rack for 5 seconds and transfer the slide rack from the first staining dish into the rinse solution within the second staining dish using forceps.

- *11.* Agitate slides within the rinse solution for 5-10 seconds and let sit at room temperature for at least 5 minutes.
- *12.* Remove the slide rack from the second staining dish. Dip into PBS-Tween20 wash buffer at room temperature and proceed with your IHC protocol.

## Water Bath Method:

- 1. Prepare Ready-to-use (RTU) solution by diluting the concentrate 1:20 in distilled water.
- 2. Fill the first staining dish with approximately 200 ml of RTU solution. Prepare a second staining dish filled with approximately 200 ml of RTU solution. The second staining dish will used later on as a rinse solution. *Staining dish size may vary, ensure sufficient quantity of RTU solution to cover the tissue slides.*
- 3. Place both staining dishes into the water bath and heat it to 97°C.
- 4. Place up to 24 formalin-fixed, paraffin-embedded tissue slides in a slide rack.
- 5. Place the slide rack into the first staining dish (ensure tissue slides are immersed in <u>pre-heated</u> RTU solution).
- 6. Incubate at 97°C for 20-40mins. The optimal incubation time should be determined by the user.
- 7. Remove both staining dishes from the water bath. Wear proper safety equipment when preforming this step.
- 8. Agitate the slide rack for 5 seconds and transfer the slide rack from the first staining dish into the rinse solution within the second staining dish using forceps.
- 9. Agitate slides within the rinse solution for 5-10 seconds and let sit at room temperature for at least 5 minutes.
- *10.* Remove the slide rack from the second staining dish. Dip into PBS-Tween20 wash buffer at room temperature and proceed with your IHC protocol.

**Alternative Protocol(s):** Accel 3-in-1 EDTA Solution is also applicable for HIER using manual pressure cookers (we recommend electric one for more consistent results) or steamers as heating sources. *The optimal conditions should be determined by the user.* 

## Precaution:

Wear gloves, proper safety equipment, and take other necessary laboratory safety precautions.

## Storage:

Store at 4°C or room temperature.

## Remarks:

For research use only.