



Polink-2 HRP Plus Rabbit DAB Detection System for Immunohistochemistry

(2-step Polymer-HRP detection system, biotin-free) **Polymer Detection System with Super Sensitivity and Specificity**

4.000	Catalog No.	D39-110	110 ml (bulk, w/o chromogen)
Storage: 4-8°C		D39-18	18 ml (with DAB, good for 180 slides)
		D39-6	6 ml (with DAB, good for 50 slides)
		☐ D39-1L	1000 ml (bulk, w/o chromogen)

Intended Use:

Polink-2 Plus HRP anti Rabbit Detection Kit is the 3rd generation of polymer detection system. It uses rabbit antibody specific enhancer to help amplify the polymer-enzyme conjugate reaction to achieve super sensitivity and specificity in immunohistochemistry staining. It produces consistent immunostaining outcomes on archival tissues and on difficult-to-work antibodies. User may need to further dilute primary antibody due to super sensitivity of Polink-2 Plus detection system. It is a biotin-free system, therefore, overcomes the non-specific staining caused by streptavidin/biotin system due to endogenous biotin. Most commonly used specimens for this system are: frozen tissue, paraffin-embedded tissue, freshly prepared lymphocytes and fixed culture cells. Suitable for manual or autostain. Staining conditions need to be optimized by user.

Polink-2 Plus HRP Detection System offers a wide choice for primary antibodies, including broad spectrum (for mouse and rabbit primary antibodies), mouse, rabbit, goat, and rat primary antibodies. Refer to **Related Product** section for details.

Kit components:

Catalog No. Product Name		Reagent 1:	Reagent 2:	Reagent3A, 3B:		
Catalog 110.	110ddct ivanic	Rabbit Antibody	Polymer HRP for Rabbit	3A : DAB Substrate (Ready-to-use)		
		Enhancer	(Ready-to-use)	3B: DAB Chromogen		
		(Ready-to-use)		(20xConcentrate)		
D39-110	Polink-2 Plus HRP DAB Rabbit Bulk kit	110ml	110ml	Not included		
D39-18	Polink-2 Plus HRP DAB Rabbit 18ml kit	18ml	18ml	30ml of DAB Reagent 3A 2ml of DAB Chromgen 3B		
				12ml of DAB Reagent 3A		
D39-6	Polink-2 Plus HRP DAB Rabbit 6ml kit	6ml	6ml	1.5ml of DAB Chromgen 3B		

Recommended Protocol:

- 1. Fixation: To ensure the quality of the staining and obtain reproducible performance, user needs to supply appropriately fixed tissue and well prepared slides.
- 2. Tissue need to be adhered to the slide tightly to avoid tissue falling off.
- 3. Paraffin embedded section must be deparaffinized with xylene and rehydrated with a graded series of ethanol before staining.
- 4. Cell smear samples should be made as much monolayer as possible to obtain satisfactory results.
- 5. Investigator needs to optimize dilution and incubation times for primary antibodies.
- Three control slides will aid the interpretation of the result: positive tissue control, reagent control (slides treated with Isotype control reagent), and negative
 control.
- 7. Staining steps: DO NOT let specimen or tissue dry from this point on.
- Serum blocking before primary antibody incubation for GBI's Polink-1, Polink-2, and Polink-2Plus is not required because all our antibody conjugates are absorbed to human serum.

Reagent	Staining Procedure	Incubation	
1. PEROXIDASE BLOCKING	a. Incubate slides in PEROXIDASE BLOCKING REAGENT (Ready-to-use 3% H ₂ O ₂ solution) for		
REAGENT. Supplied by user	10 minutes.		
	b. Rinse the slide using distilled water.		
2. HIER PRETREATMENT:	a. Heat Induced Epitope Retrieval (HIER) may be required for primary antibody suggested by		
	vendor. Please check the data sheet of primary antibody		
	b. Wash with PBS-T (PBS containing 0.05% Tween-20) 3 times for 2 minutes each time.		
3. Pre-Block (Optional)	a. Add 2 (100 µL) or more drops of 10% Normal Goat Serum (E07) to cover the tissue section and		
Supplied by user	Incubate 10 min.		
	b. Drain or blot off solution. DO NOT RINSE.		
	c. See note 8 in Recommended Protocol.		
4. PRIMARY ANTIBODY	a. Apply 2 drops (100 µL) or enough volume of PRIMARY ANTIBODY to cover the tissue section	30-60	

Supplied by user	completely. Incubate in moist chamber for 30-60 min.		
	b. Wash with PBS-T (PBS containing 0.05% Tween-20) 3 times for 2 minutes each time.		
5. Rabbit Antibody Enhancer (Ready-	a. Apply 2 drops (100 μL) or enough volume of Rabbit Antibody Enhancer to cover each section.		
to-use).	Incubate in moist chamber for 10-30 min.		
Reagent 1	(We recommend incubating the antibody enhancer up to 30mins for best sensitivity)		
	b. Wash with PBS-T (PBS containing 0.05% Tween-20) 3 times for 2 minutes each time.		
6. POLYMER-HRP for Rabbit	a. Apply 2 drops (100 µL) or enough volume of POLYMER-HRP for RABBIT to cover each	10-30	
(Ready-to-use)	section. Incubate in moist chamber for 10-30 min.		
(We recommend incubating the polymer up to 30mins for best sensitivity)			
	b. Wash with PBS-T (PBS containing 0.05% Tween-20) 3 times for 2 minutes each time.		
7. CHROMOGEN	a. Add 1 drop or 2 drops (for higher sensitivity and contrast) of Reagent 3B into 1ml of reagent 3A.		
Reagent 3A: DAB Substrate	Mix well. Protect from light and use within 5 hours.		
Reagent 3B: DAB Chromogen	b. Apply 2 drops (100 μL) or enough to completely cover tissue, of pre-mixed DAB to each section.		
	Incubate for about 5 min. Monitor the color development under microscope.		
	c. Rinse with tap water for 1-2 min.		
8. HEMATOXYLIN	a. Counterstain with 2 drops (100 ul) or enough volume of Hematoxylin to cover tissue completely	15-20	
Supplied by user	and wait about 15-20 seconds.		
	b. Rinse well under tap water for 1-2 minutes.		
	c. Put slides in PBS until show blue color (about 30-60 seconds).		
	d. Rinse well in distill or tap water		
9. Mounting medium:	Follow the manufacture data sheet procedure for mounting.		
Supplied by user	Recommended product:		
	1. GB-Mount: Cat. No. E01-18 (18ml), for AEC, AP-Red, and AP-blue.		
	2. O-Mount: Cat. No. E02-18 (18ml), for DAB		
	3. Simpo-Mount: Cat.No. E03-18 (18ml), or E03-100 (100ml), universal permanent mounting		
	medium. Can be used with or without cover slip		

Protocol Notes:

- 1. The fixation, tissue slide thickness, and primary antibody dilution and incubation time effect results significantly. Investigator needs to consider all factors and determine optimal conditions when interpret the result.
- 2. Pre-antibody blocking is optional and can be omitted if primary antibodies are diluted in buffers containing 2-10% normal goat serum.

Related Products:

Product	Catalog No.	Size	Product	Catalog No.	Size
Polink-2 Plus HRP Mouse bulk kit	D37-110	110ml	Polink-2 Plus HRP Rat-NM (No cross react	D46-110	110ml
(without chromogen)			to mouse) Bulk kit		
			(without chromogen)		
Polink-2 Plus HRP Mouse DAB kit	D37-18	18ml	Polink-2 Plus HRP Rat-NM (No cross react	D46-18	18ml
	D37-6	6ml	to mouse) DAB kit	D46-6	6ml
Polink-2 Plus HRP Goat kit Bulk Kit	D43-110	110ml	Polink-2 Plus HRP Mouse-NR (No cross	D58-110	110ml
(without chromogen)			react to Rat) Bulk kit (without chromogen)		
Polink-2 Plus HRP Goat DAB kit	D43-18	18ml	Polink-2 Plus HRP Mouse-NR (No cross	D58-18	18ml
	D43-6	6ml	react to Rat) DAB kit	D58-6	6ml
Polink-2 Plus HRP Broad Bulk Kit	D41-110	110ml	DAB Kit (2-components)	C09-12	12ml
(Without chromogen. For Ms. & Rb.)			2 tomponemo)	007 12	
					+240ml
Polink-2 Plus HRP Broad DAB Kit	D41-18	18ml	O-Mount (Organic)	E02-18	18ml
(for mouse and rabbit antibody)	D41-6	6ml			

Precautious: DAB may be carcinogenic. Please wear gloves and take other necessary precautions.

Remarks: For research use only.