


Product datasheet for CH23002-100

Netrin 1 (NTN1) Chicken Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Mouse
Host:	Chicken
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	netrin 1
Database Link:	<u>Entrez Gene 9423 Human</u>
Synonyms:	NTN1, NTN1L

Product images:




OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Dissociated cell cultures of an e13 mouse brain showing Netrin-1 (green staining) in neuronal precursor cells.




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US





Immunohistochemical colocalization of netrin-1 and endothelial markers. Paraformaldehvdefixed frozen tissue sections were stained with anti-netrin-1 antibody (red stain) and anti-Von Willebrand Factor (VWF) antibody (green stain) as described in MATERIALS AND METHODS. Nuclei were stained blue with DAPI. In sham-operated kidney (B-D), netrin-1 staining was mainly in the interstitium and colocalized with VWF. Similarly, netrin-1 and VWF colocalized in the kidney section from 72 h of reperfusion (E and F). No staining was seen in the primary antibody negative control (A). Colocalization of netrin-1 (red) and GFP (green) in kidney from Tie2-GFP mice was also observed under normal conditions (G and H). doi:10.1152/ajprenal.00508.2007.


Immunohistochemical localization of netrin-1 in renal cell carcinoma (RCC) tissues. A. Secondary antibody control showing no staining. B. Normal adjacent tissues do not show any staining for netrin-1. C. Stage I RCC shows staining for netrin-1. D. Stage II RCC shows staining for netrin-1. E&F. doi: 10.3109/1354750X.2010.541564



Dissociated cell cultures of an e13 mouse brain showing Netrin-1 (green staining) in neuronal precursor cells.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US





Immunohistochemical colocalization of netrin-1 and endothelial markers. Paraformaldehydefixed frozen tissue sections were stained with anti-netrin-1 antibody (red stain) and anti-Von Willebrand Factor (VWF) antibody (green stain) as described in MATERIALS AND METHODS. Nuclei were stained blue with DAPI. In sham-operated kidney (B-D), netrin-1 staining was mainly in the interstitium and colocalized with VWF. Similarly, netrin-1 and VWF colocalized in the kidney section from 72 h of reperfusion (E and F). No staining was seen in the primary antibody negative control (A). Colocalization of netrin-1 (red) and GFP (green) in kidney from Tie2-GFP mice was also observed under normal conditions (G and H). doi:10.1152/ajprenal.00508.2007.

Immunohistochemical localization of netrin-1 in renal cell carcinoma (RCC) tissues. A. Secondary antibody control showing no staining. B. Normal adjacent tissues do not show any staining for netrin-1. C. Stage I RCC shows staining for netrin-1. D. Stage II RCC shows staining for netrin-1. E&F. doi: 10.3109/1354750X.2010.541564

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US