


Product datasheet for **EB200001**

Lenti RapidTiter Sticks

Product data:

Product Type:	Lenti Titer Kits
Description:	Lateral flow testing sticks used for the semi-quantitative detection of the lentiviral p24 protein, 5 tests
Symbol:	HIV-1 p24
Background:	The gag protein p24 (MW 24kD) is a major structural protein of the HIV capsid, which is required for virus particle assembly. There are approximately 2000 p24 molecules per virus particle, or at a molecule weight of 24 kDa, about 104 virus particles per picogram of p24. HIV-1 p24 measurement is commonly used as an indicator of HIV-1 infection and viral load. Quantitative detection of the HIV-1 p24 has become a standard method for measuring lentiviral titer.
Storage:	Upon receipt, store at room temperature
Stability:	The kit should not be used beyond the expiration date. Leave Lenti RapidTiter Sticks sealed in its foil pouch until just before use. Do not use if pouch is damaged or open. Once opened, the Lenti RapidTiter Stick should be used immediately.

Product images:



[View online »](#)