

PreBlock Solution
(10% Normal Goat Serum)

Storage: 2-8°C

Catalog No:	E07-18	18 ml
	E07-100	100ml
	E07-500	500ml



Intended Use:

PreBlock Solution is 10% normal goat serum solution. Normal goat serum is used to block non-specific staining when a goat second antibody is used during immunohistochemical staining procedure.

Kit Components:

Catalog No.	Size
E07-18	18 ml (ready to use)
E07-100	100ml (ready to use)
E07-500	500ml (ready to use)

Solution contains 0.05% NaN₃.

Recommended Protocol:

1. Prior to application of primary antibody apply PreBlock Solution on tissue slides and incubate the tissue section for 10 minutes at room temperature
2. Blot off PreBlock Solution (do not rinse) and add primary antibody

Precaution:

Wear gloves and take other necessary laboratory safety procedure.

Storage:

Store at 2-8°C.

Remarks:

For research use only.