

Product datasheet for RC215036L3V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ABCG1 (NM_207174) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: ABCG1 (NM_207174) Human Tagged ORF Clone Lentiviral Particle

Symbol: ABCG1

Synonyms: ABC8; WHITE1

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

 Tag:
 Myc-DDK

 ACCN:
 NM_207174

ORF Size: 2031 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC215036).

OTI Disclaimer:

Sequence:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 207174.1, NP 997057.1

 RefSeq Size:
 2983 bp

 RefSeq ORF:
 2034 bp

 Locus ID:
 9619

 UniProt ID:
 P45844

 Cytogenetics:
 21q22.3

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: ABC transporters





ORÏGENE

MW: 75.2 kDa

Gene Summary: The protein encoded by this gene is a member of the superfamily of ATP-binding cassette

(ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the White subfamily. It is involved in macrophage cholesterol and phospholipids transport, and may regulate cellular lipid

homeostasis in other cell types. Six alternative splice variants have been identified. [provided

by RefSeq, Jul 2008]