

Product datasheet for RG217916

TARP (NM_001003806) Human Tagged ORF Clone

Product data:

OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Expression Plasmids
21	
Product Name:	TARP (NM_001003806) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TARP
Synonyms:	CD3G; TCRG; TCRGC1; TCRGC2; TCRGV
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	<pre>>RG217916 representing NM_001003806 Red=Cloning site Blue=ORF Green=Tags(s)</pre>
	TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCC <mark>GCGATCGC</mark> C
	ATGAAGACTAACGACACATACATGAAATTTAGCTGGTTAACGGTGCCAGAAAAGTCACTGGACAAAGAAC ACAGATGTATCGTCAGACATGAGAATAATAAAAACGGAGTTGATCAAGAAATTATCTTTCCTCCAATAAA GACAGATGTCATCACAATGGATCCCAAAGACAATTGTTCAAAAGATGCAAATGATACACTACTGCTGCAG CTCACAAACACCTCTGCATATTACATGTACCTCCTCCTGCTCCTCAAGAGTGTGGTCTATTTTGCCATCA TCACCTGCTGTCTGCTTAGAAGAACGGCTTTCTGCTGCAATGGAGAGAAATCA
	ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA
Protein Sequence:	>RG217916 representing NM_001003806 <mark>Red=</mark> Cloning site Green=Tags(s)
	MKTNDTYMKFSWLTVPEKSLDKEHRCIVRHENNKNGVDQEIIFPPIKTDVITMDPKDNCSKDANDTLLLQ LTNTSAYYMYLLLLLKSVVYFAIITCCLLRRTAFCCNGEKS
	TRTRPLE - GFP Tag - V
Restriction Sites:	Sgfl-Mlul




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Cloning Scheme:



Plasmid Map:



ACCN:	
ORF Size:	
OTI Disclaimer:	

NM_001003806

333 bp

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Service TARP (NM_001003806) Human Tagged ORF Clone – RG217916	
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	 Centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM 001003806.1, NP 001003806.1</u>
RefSeq Size:	1027 bp
RefSeq ORF:	336 bp
Locus ID:	445347
Cytogenetics:	7p14.1
Protein Families:	Transmembrane
Gene Summary:	In some non-lymphoid tissues, the unrearranged T cell receptor gamma (TRG@) locus is expressed. The resulting transcript contains a subset of the TRG@ gene segments and is shorter than TRG@ transcripts expressed in lymphoid tissues. This RefSeq record represents the unrearranged TRG@ locus transcript; the complete TRG@ locus is represented by the genomic RefSeq NG_001336. The transcript represented by this RefSeq has two open reading frames (ORFs) that encode different proteins. The downstream ORF is in the same frame as TRG@ and its protein product is similar to TRG@ proteins. The upstream ORF uses a different

reading frame and encodes a novel protein. [provided by RefSeq, Jul 2008]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US