

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Product datasheet for SC336634

## Dystrophia myotonica protein kinase (DMPK) (NM\_001288766) Human Untagged Clone


## **Product data:**

Expression Plasmids
Dystrophia myotonica protein kinase (DMPK) (NM_001288766) Human Untagged Clone
Tag Free
DMPK
DM; DM1; DM1PK; DMK; MDPK; MT-PK
Neomycin
pCMV6-Entry (PS100001)
Kanamycin (25 ug/mL)



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Dystrophia myotonica protein kinase (DMPK) (NM_001288766) Human Untagged Clone – SC336634
Fully Sequenced	<b>ORF:</b> >SC336634 representing NM_001288766.
i any bequenced	Blue=Insert sequence Red=Cloning site Green=Tag(s)
	GCTCGTTTAGTGAACCGTCAGAATTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
	GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
	ATGTCAGCCGAGGTGCGGCTGAGGCGGCTCCAGCAGCTGGTGTTGGACCCGGGCTTCCTGGGGCTGGAG
	CCCCTGCTCGACCTTCTCCTGGGCGTCCACCAGGAGCTGGGCGCCTCCGAACTGGCCCAGGACAAGTAC
	GTGGCCGACTTCTTGCAGTGGGCGGAGCCCATCGTGGTGAGGCTTAAGGAGGTCCGACTGCAGAGGGAC
	GACTTCGAGATTCTGAAGGTGATCGGACGCGGGGCGTTCAGCGAGGTAGCGGTAGTGAAGATGAAGCAG
	ACGGGCCAGGTGTATGCCATGAAGATCATGAACAAGTGGGACATGCTGAAGAGGGGCGAGGTGTCGTGC
	TTCCGTGAGGAGAGGGACGTGTTGGTGAATGGGGACCGGCGGTGGATCACGCAGCTGCACTTCGCCTTC
	CAGGATGAGAACTACCTGTACCTGGTCATGGAGTATTACGTGGGCGGGGACCTGCTGACACTGCTGAGC
	AAGTTTGGGGAGCGGATTCCGGCCGAGATGGCGCGCGTTCTACCTGGCGGAGATTGTCATGGCCATAGAC
	TCGGTGCACCGGCTTGGCTACGTGCACAGGGACATCAAACCCGACAACATCCTGCTGGACCGCTGTGGC
	CACATCCGCCTGGCCGACTTCGGCTCTTGCCTCAAGCTGCGGGCAGATGGAACGGTGCGGTCGCTGGTG
	GCTGTGGGCACCCCAGACTACCTGTCCCCCGAGATCCTGCAGGCTGTGGGCGGTGGGCCTGGGACAGGC
	AGCTACGGGCCCGAGTGTGACTGGTGGGCGCTGGGTGTATTCGCCTATGAAATGTTCTATGGGCAGACG
	CCCTTCTACGCGGATTCCACGGCGGAGACCTATGGCAAGATCGTCCACTACAAGGAGCACCTCTCTCT
	CCGCTGGTGGACGAAGGGGTCCCTGAGGAGGCTCGAGACTTCATTCA
	ACACGGCTGGGCCGGGGTGGAGCAGGCGACTTCCGGACACATCCCTTCTTCTTTGGCCTCGACTGGGAT
	GGTCTCCGGGACAGCGTGCCCCCTTTACACCGGATTTCGAAGGTGCCACCGACACATGCAACTTCGAC
	TTGGTGGAGGACGGGCTCACTGCCATGGAGACACTGTCGGACATTCGGGAAGGTGCGCCGCTAGGGGGTC
	CACCTGCCTTTTGTGGGCTACTCCTACTCCTGCATGGCCCTCAGGGACAGTGAGGTCCCAGGCCCCACA
	CCCATGGAACTGGAGGCCGAGCAGCTGCTTGAGCCACACGTGCAAGCGCCCAGCCTGGAGCCCTCGGTG
	TCCCCACAGGATGAAACAGCTGAAGTGGCAGTTCCAGCGGCTGTCCCTGCGGCAGAGGCTGAGGCCGAG
	GTGACGCTGCGGGAGCTCCAGGAAGCCCTGGAGGAGGAGGTGCTCACCCGGCAGAGCCTGAGCCGGGAG
	ATGGAGGCCATCCGCACGGACAACCAGAACTTCGCCAGTCAACTACGCGAGGCAGAGGCTCGGAACCGG
	GACCTAGAGGCACACGTCCGGCAGTTGCAGGAGCGGATGGAGTTGCTGCAGGCAG
	CCCTAG
	ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
	TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
Restriction Sites:	Sgfl-Mlul
Plasmid Map:	
	fl ori



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	strophia myotonica protein kinase (DMPK) (NM_001288766) Human Untagged Clone – 336634
ACCN:	NM_001288766
Insert Size:	1593 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Meth	<ul> <li>od: 1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ul>
RefSeq:	<u>NM 001288766.1</u>
RefSeq Size:	2722 bp
RefSeq ORF:	1593 bp
Locus ID:	1760
Cytogenetics:	19q13.32
Protein Families:	Druggable Genome, Protein Kinase
MW:	59.4 kDa

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Dystrophia myotonica protein kinase (DMPK) (NM_001288766) Human Untagged Clone – SC336634
Gene Summary:	The protein encoded by this gene is a serine-threonine kinase that is closely related to other kinases that interact with members of the Rho family of small GTPases. Substrates for this enzyme include myogenin, the beta-subunit of the L-type calcium channels, and phospholemman. The 3' untranslated region of this gene contains 5-38 copies of a CTG trinucleotide repeat. Expansion of this unstable motif to 50-5,000 copies causes myotonic dystrophy type I, which increases in severity with increasing repeat element copy number. Repeat expansion is associated with condensation of local chromatin structure that disrupts the expression of genes in this region. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Jul 2016] Transcript Variant: This variant (7) has multiple differences in the presence and absence of exons at its 5' end, compared to variant 1. These differences produce a distinct 5' UTR, cause translation initiation at an alternative start codon, loss of an in-frame portion of the 3' coding region, and a frameshift in the 3' coding region, compared to variant 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US